


# NASA TECH BRIEF



NASA Tech Briefs are issued by the Technology Utilization Division to summarize specific technical innovations derived from the space program. Copies are available to the public from the Clearinghouse for Federal Scientific and Technical Information, Springfield, Virginia, 22151.

## Rotating Holder Permits Accurate Grinding of Metallurgical Microsamples



**The problem:** To accurately grind or polish a level, flat surface on metallurgical microsamples. Previous methods held the sample in a fixed position and moved it in a straight line across a rotating grinding wheel. This resulted in finished surfaces with a slight slope (not exactly perpendicular to the center axis of the mounted specimen).

**The solution:** A fixture that rotates the sample approximately 180° with each pass across the grinding wheel. This rotation minimizes any sloping of

the sample surface by irregularities in the grinding wheel.

**How it's done:** The sample is mounted rigidly in the holder and the holder is soldered to a rod that mates by a sliding fit to a mounting frame. The rod is free to rotate in the mounting frame and, as the sample passes across the rotating grinding wheel, the holder and sample turn approximately 180° in the mounting frame. The direction of grinding on one side of the wheel is reversed as the sample moves to the opposite side of the wheel. The weight of

(continued overleaf)

holder and rod applies the necessary pressure for grinding.

**Notes:**

1. This device could be modified to grind a number of samples simultaneously.

2. Inquiries concerning this innovation may be directed to:

Technology Utilization Officer  
Lewis Research Center  
21000 Brookpark Road  
Cleveland, Ohio, 44135  
Reference: B65-10262

**Patent status:** NASA encourages commercial use of this innovation. No patent action is contemplated by NASA.

Source: Donald L. Cramer  
(Lewis-131)